



**Planmed**

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Planmed	Nuance Manager 3	3.5	May 2013
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query	
		Broad Worklist Query	
		PPS Exception Management	
		Modality Group Case	
Mammography Image	Acquisition Modality	Partial View	
Consistent Time	Time Client	None	
Patient Information Reconciliation	Acquisition Modality	None	
Internet address for vendor's IHE information: <a href="http://www.planmed.com/ihe">http://www.planmed.com/ihe</a>			
Links to Standards Conformance Statements for the Implementation			
<b>HL7</b>	not applicable		
<b>DICOM</b>	<a href="http://www.planmed.com/dicom">http://www.planmed.com/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	